

3400 & 3409

Wall Mounted, Telescopic Screen



CUSTOMER INSTRUCTIONS

PLEASE READ AND FAMILIARIZE YOURSELF WITH ALL INSTRUCTIONS BEFORE USING THIS PRODUCT. If you have trouble understanding these instructions contact your dealer or Winco customer support, (800) 237-3377, e-mail: customerservice@wincomfg.com before attempting to use this product; otherwise injury may occur.

Winco assumes no responsibility for damage or injury caused by the improper assembly, installation or use of these products.



1. READ AND FOLLOW ALL DIRECTIONS.
2. NEVER hang equipment or personal objects from Curtain Rod; Screen could fall and INJURY MAY OCCUR.
3. NEVER place screen near open flame or over a radiator or heat register.
4. NEVER place screen in front of a door or exit sign; In case of emergency screen could block or hinder evacuation.
5. SUFFOCATION DANGER to children; Screen could block nose and mouth and prevent breathing. VINYL SCREEN IS NOT A TOY.
6. DO NOT allow children to play with, or climb on screen; Screen could fall and SERIOUS INJURY MAY OCCUR.
7. Periodically, recheck tightness of mounting screws.
8. ALWAYS mount at least two mounting screws to a wall stud, "dead wood" or other solid wall structure; DO NOT attach to drywall; Professional installation is recommended.
9. DO NOT use screen in or with ANY type of vehicle or trailer.
10. Immediately remove from service; Any screen with broken or missing parts, torn vinyl, or other visible damage.
11. USE ONLY WINCO AUTHORIZED REPLACEMENT PARTS.
12. SAVE THESE INSTRUCTIONS for future reference and training.

Tools needed:


- Phillips-Head Screwdriver
- Scissors or box cutter
- Special tools for mounting to concrete wall (ex. Drill, masonry bit, etc...) it is recommended that someone with experience using these types of tools perform this type of installation.
- - If using powered tools always wear safety goggles –

PLEASE READ AND FAMILIARIZE YOURSELF WITH ALL INSTRUCTIONS BEFORE PROCEEDING WITH ASSEMBLY

MOUNTING BRACKET TO WALL

1. You must first choose a sturdy, unobstructed (clear of ceiling, fans, electrical wiring, etc...) mounting location approximately 80" from floor to top of Mounting Bracket.
IMPORTANT: Mounting to drywall only is NOT recommended.
2. If mounting to a wall with a solid wood backer-board; Use all four (4) Phillips screws provided; DO NOT OVERTIGHTEN.
3. If mounting to a cement or concrete-block wall use toggle bolts or appropriate concrete anchors – capable of at least 30 lbs. holding strength.
4. Install (4) screws – through holes in bracket and into wall – tighten securely & snap on plastic cover.

SEE PHOTOS – NEXT PAGE



INSTALL ROD & CURTAIN

5. Place short end of Rod Assmly through both holes on Wall Bracket; Nylon washer should make contact with top of wall bracket to give the 340 screen smooth operation when rotating.
6. "Unsnap" curtain hook latches and attach curtain hooks through holes in Vinyl Panel, then around Rod Assembly.
7. Secure hooks to Rod Assembly by "snapping" the hook latches.


PANEL REPLACEMENT

8. To replace Vinyl Panels, remove the curtaining hooks from the vinyl panel by "unsnapping" hook latch.
9. Replace old vinyl panel with a new vinyl panel and re-install curtain hooks.

CARE & MAINTENANCE OF YOUR SCREEN

10. PERIODICALLY CHECK ALL MOUNTING HARDWARE FOR TIGHTNESS

Wrinkles in Vinyl – Storage and shipping causes wrinkles of the vinyl of your Winco Screen; This is normal. Please allow 1-2 weeks at room temperature (approx. 70+ degrees) to allow for vinyl to return to original appearance.



Visit our website www.wincomfg.com/product-documents/parts-lists for a complete parts list catalog of your chair or call Customer Service @ 1-800-237-3377 to request it by mail.

GENERAL CLEANING:

Winco furnishings are built to provide durable reliable service when properly cared for. In general all of our products should be 1) Cleaned 2) Rinsed 3) Allowed to air dry. It is important to note that disinfecting a product is not the same as cleaning the product. Disinfectants alone will not provide adequate cleaning, they do not have grease or oil cutting properties to remove grime or hair and skin oils. Over time, the appearance and feel and performance of your upholstery may diminish if not cleaned properly

CLEANING:

Remove grime, hair and body oils; Your furnishings should be cleaned with mild soap & water solution and a damp cloth on a regular basis. Avoid harsh detergents or chemicals that could damage the finish of your product. If the furnishing is disinfected with chemicals or bleach, it **MUST** be wiped off using only clean water on a damp cloth and then be allowed to air dry. Failure to rinse with clean water can result in a build-up of residues over time that can damage your product.





LIMITED WARRANTY

• 5516 SW 1st Lane • Ocala, FL 34474
 • Phone (352) 854-2929 • (800) 237-3377
 • Fax (352) 854-9544
 • www.wincomfg.com
 • customerservice@wincomfg.com

TERMS

Winco Mfg., LLC. warrants this product to be free of manufacturer's defects in material and workmanship, provided that the product is used according to normal operating conditions and proper maintenance intended by manufacturer. This warranty is available only to the original retail purchaser, is non-transferable and commences on the date of retail sale; proof of purchase required.

Limited One Year Product Warranty: Winco Mfg., LLC warrants the complete product for one (1) year. At Winco Mfg., LLC. sole discretion, it may repair or replace any components freight free that are found to be defective during the first year. *Winco Mfg., LLC. shall not be liable for any labor, or any other costs incurred as a result of or in conjunction with a warranty claim.*

Limited Three Year Warranty: Winco Mfg., LLC warrants recliner mechanisms, electrical components, vinyl panels and Vinyl on upholstered parts [from cracking or delaminating] pressurized gas springs, and casters, when new, for a period of three (3) years. At Winco Mfg., LLC. 's Sole discretion, it may repair or replace components found to be defective. *Winco Mfg., LLC. shall not be liable for any labor, shipping or any other costs incurred as a result of or in conjunction with a warranty claim.*

Limited Lifetime Frame Warranty: Winco Mfg., LLC warrants steel and aluminum base frames for all products, when new, for the lifetime of the original purchaser. This limited warranty does not apply to paint/finish or any components attached to the frame such as; upholstery, foam, casters, mechanisms or related parts that are covered under above warranties. *Winco Mfg., LLC shall not be liable for any labor, shipping or any other costs incurred as a result of or in conjunction with a warranty claim.*

Design, appearance, parts, & construction are subject to change without notice

The purchaser's exclusive remedy under this warranty shall be limited to such repair or replacement of defective components at Winco Mfg., LLC sole discretion. For warranty service, contact the dealer that the product was originally purchased from or Winco Mfg., LLC directly.

EXCLUSIONS

There are no other warranties, conditions, representations or guarantees, express or implied, made or intended by Winco Mfg., LLC and all other warranties, conditions, representations or guarantees including any warranties, conditions, representations or guarantees under any Sale of Good Act or Like legislation or statute is hereby expressly excluded. Any and all other implied warranties shall not extend beyond the duration of the express warranty. Liability for incidental or consequential damages is excluded to the extent permitted by law. Some states do not allow the exclusion or limitation so the above limitation or exclusion may not apply to you. This warranty gives you specific legal rights. In addition you may also have other rights, which vary from state to state.

GENERAL PROVISIONS

No warranty herein contained or set out shall apply when damage or repair is caused by any of the following:

1. Damage in transit.
2. Accident, alteration, abuse or misuse of product i.e. exceeding weight capacities, applying inordinate pressure to footrest/leg rest, use in Mental Health facilities like institutions or, any unintended use of the product or use in unintended environments (i.e. outdoor, showers, MRI rooms, etc..)
3. Fire, water damage, theft, war, riot, hostility, acts of God.

Examples contained in this list are not to be construed as all-inclusive.

