

SUGGESTED PROTOCOL FOR CHARGING YOUR TMM BATTERY:

1. Ensure the cord for your **TMM Battery Charger** is fully inserted into the battery charger inlet.
2. Plug your **TMM Battery Charger** directly into an outlet. Do **NOT** use an extension cord or power strip. Extension cords and/or power strips can cause a drop in current that may hinder the charger's ability to charge the battery efficiently.
3. Slide the **TMM Battery** downward onto the **TMM Battery Charger** making sure the battery is properly seated onto the charging port.
4. Check to see if either an orange light (indicating the battery is "charging") or a green light (indicating the battery is "fully charged") appears.
5. Charge your **TMM Battery** for eight (8) hours or until the Green light indicator appears to ensure an optimum charge is achieved.
6. Rotate batteries from charger to chair every two (2) to three (3) days to maintain an optimum charge on the batteries.
7. When your TMM Stretcher-Chair is not in use for extended periods of time (four (4) days or more) unplug the battery from the Stretcher-Chair to prevent battery drainage.
8. Replace your **TMM Battery** every two (2) years (recommended) or as required.