



SUGGESTED MAINTENANCE SCHEDULE

TREATMENT AND CLINICAL RECLINERS

DAILY:

Clean according to facility requirements, and RINSE chair with clean water according to User Instruction Manual (available online and specific to model).

WEEKLY:

Place chair in Trendelenburg position and ensure proper operation
Visual inspection for any loose, torn, cracked vinyl or hardware components.

MONTHLY:

Ensure Heat & Massage is working. (No lights on switch are flashing)
Ensure there is no cut, torn or cracked vinyl.
Ensure chair mechanism operates freely, full recline to full upright positions.
Ensure that all casters roll freely and lock correctly.
Extensive overall cleaning recommended (wipe down frame, casters underside etc.)

QUARTERLY:

Inspect all electrical components, ensure:

- Cables are routed properly (No kinks, binding or rubbing)
- There are no breaks, cracks, cuts or exposed wires
- Inspect Line cord (power) where it enters the control box
- With Swing-away arm models, Insure Arm latching components are tight and aligned, and functioning properly.

Ensure swing-away arm hinges are tight and aligned

Inspect any cables (Backrest release, Trendelenburg release, etc.) for any signs of fraying or binding.

Inspect side tables for proper operation

SEMI ANNUAL:

Lubricate all pivot points on mechanism with facility approved lubricant.

Check integrity and tightness of all hardware (Nuts, bolts, screws.)

ANNUAL

Inspect frame for integrity, touch up scratches and inspect for corrosion.

CONTACT US

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Insidesales@WincoMfg.com

MS149908 REV A