




# Vinyl Selections




02 ROYAL BLUE



03 TAUPE



04 MOSS GREEN




06 HUNTER GREEN



14 BURGUNDY



16 MAUVE



17 BLUERIDGE




07 GREY



60 PYRITE



61 GYPSUM



62 ICE BLUE




63 CYPRESS



34 AUBERGINE



32 GROTTA



33 EUCALYPTUS



37 NAVY



41 BRIAR



42 ROSEBUD



43 EVERGREEN



44 NIGHTSHADE

**Vinyl colors displayed here may vary slightly from actual color.  
Actual vinyl samples should always be viewed before placing an order.**

*Additional optional vinyl colors available. Please contact Winco for information.*

- Upholstery material is heavy gauge, expanded vinyl that is easily cleaned.
- The last two digits of product code indicates upholstery color. Example 655N03 is a Drop Arm Care Cliner in Taupe.
- All optional vinyls still contain the same contract vinyl specifications as our standard vinyls.
- All foam used in the upholstery products is medium density polyurethane, unless otherwise noted. It is certified to meet Motor Vehicle Safety Standard #302 and California Tech. Bulletin #117
- TB-133 compliance option available at a nominal upcharge.
- Vinyl meets the following specifications: Federal Aviation Specification 25.853/Federal Specification CCC-A-680A/Boston Fire Department Code.



09 WHITE *Table & Stools Only*



08 BLACK *Table & Stools Only*

These colors are standard on model #556/ 656/690