

MODEL OPTIONS

STRETCHER-CHAIRS

REFERENCE SHEET

OPTIONS are selected at the time of purchase and are mounted/welded onto your Stretcher-Chair at our factory in Ocala, FL.

Letters are used to represent each of the options available (as shown on the right). The letters for the options you choose will appear in reverse alphabetical order after the chair model you select (For Example: TMM4PLUS-XWTRFEB)

PLEASE NOTE: all of the options shown may not be available for every model and/or in combination with one another.

Contact your TransMotion Representative to discuss your specific needs and to learn more about which options will provide you with the best solution.

			
X: SEAT TILT	W: WIDE	T: TRIPLE COLUMN	S: SHORT BASE
			
R: UPPER RAILS	Q: SCALE	L: LOW HEIGHT	G: ENG BACK
			
F: FOLDING FOOTREST	E: EXPORT PLUG	B: BATTERY	A: 120V A/C

ONE PATIENT, ONE SURFACE™

 **TransMotion**
by Winco

TM-147141 REV B