

CONVALESCENT RECLINER

5251 (WITH TRAY) / 5261 (WITHOUT TRAY)

5 SERIES



5251 Shown

Practicality, Durability and Comfort...

- Built-in Headrest
- Side storage pocket
- Recessed hand-grips
- Push handle for safer mobility
- 3 Locking positions: upright, TV or recline
- Model 5251: non-porous tray stows at side of chair when not in use

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



SPECIFICATIONS

CONVALESCENT RECLINER	5251	5261
DESCRIPTIONS	TRAY INCLUDED	NO TRAY
Seat Height	19.5" (49.5 cm)	
Seat Width	19" (48.3 cm)	
Seat Depth	19.5" (49.5 cm)	
Overall Height	45" (114.3 cm)	
Overall Width	26.5" (67.3 cm)	
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	33" (83.8 cm)	
Armrest Height	27" (68.6 cm)	
Overall Length (Back Upright)	39" (99.1 cm)	
Overall Length (Back Reclined)	64" (162.6 cm)	
Weight	93 lbs (42.2 kg)	88 lbs (39.9 kg)
Weight Capacity	275lbs (124.7 kg)	

SUGGESTED OPTIONS & ACCESSORIES



IV POLE
Left, Right, or
Dual



**FLIP-UP
FOOTRESTS**

OVER 50
DESIGNER
COLORS
AVAILABLE

Additional options and accessories available.

Please visit our website to view a full selection of options or contact your local representative.

WARRANTY Contact Winco or visit our website for warranty details.

COLORS Visit our website or contact Winco for available color selections and/or swatches

NOTE Winco reserves the right to change specifications without notice.

CERTIFICATIONS 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com