

# DROP-ARM CONVALESCENT RECLINER

5271 (WITH TRAY) / 5281 (WITHOUT TRAY)

5 SERIES



5281 Shown

## Practicality, Durability and Comfort Combined with Easy-Access...

- Dual Drop-Arms with One-Touch Release
- Built-in Headrest
- Padded Push Handle
- 3 Locking positions: upright, TV or recline
- Model 5271: non-porous tray stows at side of chair when not in use


ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# SPECIFICATIONS

DROP-ARM CONVALESCENT RECLINER	5271	5281
DESCRIPTIONS	TRAY INCLUDED	NO TRAY
Seat Height	21.5" (54.6 cm)	
Seat Width	19" (48.3 cm)	
Seat Depth	21.75" (55.2 cm)	
Overall Height	46" (116.8 cm)	
Overall Width	26" (66 cm)	
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	32" (81.3 cm)	
Overall Length (Back Upright)	41" (111.8 cm)	
Overall Length (Back Reclined)	65" (165.1 cm)	
Weight	93 lbs (42.2 kg)	88 lbs (39.9 kg)
Weight Capacity	275lbs (124.7 kg)	

## SUGGESTED OPTIONS & ACCESSORIES



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**