

# LIFE CARE CLINER

5 SERIES

5851 (WITH TRAY) / 5861 (WITHOUT TRAY)



## Value, Durability and Comfort...

- Molded side-panels
- Push handle for easy maneuverability
- 3 Locking positions: upright, TV or recline
- Padded armrests with scuff guards
- Recessed handgrips provide additional assistance for ingress and egress
- Model 5851: Side storage for lap tray


ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



## SPECIFICATIONS

| LIFE CARE CLINER               | 5851              | 5861             |
|--------------------------------|-------------------|------------------|
| DESCRIPTIONS                   | TRAY INCLUDED     | NO TRAY          |
| Seat Height                    | 20" (50.8 cm)     |                  |
| Seat Width                     | 19" (48.3 cm)     |                  |
| Seat Depth                     | 20" (50.8 cm)     |                  |
| Overall Height                 | 46.5" (118.1 cm)  |                  |
| Overall Width                  | 26.5" (67.3 cm)   |                  |
| Range of Seat Angle            | 8° - 32°          |                  |
| Range of Back Angle            | 62° - 38°         |                  |
| Back Height above Seat         | 32" (81.3 cm)     |                  |
| Overall Length (Back Upright)  | 41" (104.1 cm)    |                  |
| Overall Length (Back Reclined) | 63" (160 cm)      |                  |
| Weight                         | 85 lbs (38.6 kg)  | 80 lbs (36.3 kg) |
| Weight Capacity                | 275lbs (124.7 kg) |                  |

## SUGGESTED OPTIONS & ACCESSORIES



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**