



PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1	A	404945	SEAT PLATE W/LEVER FOR TASK STOOL	1
1A		566517	SCREW-MACH 1/4-20 X 1 PAN PHILLIP M/S ZINC	4
1B		566303	WASHER-1/4 SAE FLAT WASHER ZINC	4
2		405000	22"BLK STOOL BASE 50MM (430/435/438)	1
2A		405200	GAS CYLINDER 140MM STROKE 430/435	1
3A	B	405100	2" 50mm CASTER (430/435/438)	5
3B	C	405120	Caster,High Resistance Braking, Nylon Twin Wheel	5
3C		405120K	Caster,High Resistance Braking, Nylon Twin Wheel KIT	1

NOTES

Note	Description
A	For 435's WITH an "A" in the serial # ONLY (EXAMPLE: 435A555555)
B	4350 Stool Models.
C	4358 Stool Models



PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1	C	U0348XX	4358 ONLY UPHOLSTERED SEAT	1
1	A	U0341XX	435/438/445 UPHOLSTERED SEAT	1
2	A,B	U0240XX	435A/438A UPHOLSTERED BACK	1

NOTES

Note	Description
A	XX- Use color code for specific part number i.e. 02 - Royal Blue For 435's WITH an "A" in the serial # ONLY (EXAMPLE: 435A555555).
B	Upholstery ONLY not available seperately. Part ships assembled.
C	For 4358 stool models ONLY.
OPTION	For TB133 Add TB after Part Number i.e. U010702 TB
	Use 00 in place of options codes NOT applicable i.e. U100602 00 For Uph Seat without TB

REV	Effective Date	CCR Number	Description of Change
01	04/12/2017	17-049	Add white vinyl backing to all 4358 stools