



Clinical staff, surgeons, and patients can all benefit when a Stretcher-Chair is used within multiple healthcare settings.

OPERATIONAL AND STAFF EFFICIENCY:

- Convertible equipment streamlines the process, saving time and space, by allowing patients to go from admission, to pre-op, to procedure, to discharge on a single device

IMPROVED SAFETY:

- Limiting the number of times a patient must be lifted, re-positioned, or transferred reduces the risk of staff injuries
- Enabling patients to remain on one surface throughout an out-patient stay lowers the potential for fall related injuries

ENHANCED CUSTOMER SATISFACTION:

- A range of positioning options eases transport of patients who are experiencing discomfort prior to or preceding their procedure
- Enabling patients to remain on one surface throughout an out-patient stay lowers the potential for fall related injuries
- TMM Series Stretcher-Chairs are equipped with 3.75" Memory Foam cushions, creating a more comfortable experience

TMM4 ADVANTAGE:

- Precise positioning at the push of a button
- 500lb weight capacity and 24" seat width accommodates a wider range of patients
- Full radiolucent back allows for in-unit X-ray procedures
- Self-storing push bars fold down for an unobstructed view of your patient
- Optional integrated scale allows for weight measurements to be taken while the patient is seated
- Brake and steer caster system for easy mobility
- Configurable to meet your specific needs - select from an assortment of optional enhancements and accessories