

Stretcher-Chairs streamline the treatment process and surgical flow by replacing multiple pieces of equipment with a single device.

PRE-OP



TRANSPORT



PROCEDURE



RECOVERY



## TMM5 ADVANTAGE:

- Accurate temporal positioning for effective surgical access
- Reduced need to transfer patients limits the risk of injury to patients and staff
- Easily position the patient with the push of a button

## KEY FEATURES:

- Low profile, dual articulating headrest tilts from 30° to 60° as well as raised or lowered
- Adjustable height (25" - 33" or 25" - 40") provides ergonomic access to patient
- 500lb weight capacity
- 24" or 28" seat widths to accommodate a wider range of patients
- 3.75" Memory Foam cushions help reduce pressure injuries and improve patient comfort
- Self-storing push bars fold down for an unobstructed view of your patient
- Configurable to meet your specific needs - select from an assortment of optional enhancements and accessories. For Ophthalmic procedures, we suggest the TMM5-XTB (shown above).

